

How to get technical support?

Troubleshooting Clarity Lite

If you encounter any problems while using Clarity Lite, please do not hesitate to contact technical support for help via email at support@dataapex.com. To diagnose your problem, DataApex will need files containing information about the settings of your Clarity Lite station and a record of the last performed actions (a complete list of requested files is provided below).

Before you contact technical support

- Consult online help
Pressing the **F1** key in the dialog where you encountered the problem will invoke **context help**.
- Check our website for **available program updates** to your version – only the latest version is supported
Use the *Help – Check for Updates* command in the Clarity Lite window (internet connection necessary) or view the *Downloads – Full versions* section at www.dataapex.com.

Collecting information for technical support

When contacting [DataApex technical support](mailto:support@dataapex.com), please provide the **description of your problem** and send us the following files (as sorted by priority from the highest):

1. Files from the Clarity Lite installation folder (C:\Clarity_Lite\Cfg)

- SystemInfo.txt, Clarity.cfg and all *.dsk files

2. Files from the Log folder (C:\Clarity_Lite\Cfg\Audit_Trails)

- YYYY_MM_DD.audit files from the day when the issue occurred + from a few days before and after
- YYYY_MM_DD.audit.bak (might not be present)
- Global audit trail from the subfolder \Global (YYYY_MM_DD_global.audit)

3. Everything from the subfolder C:\Clarity_Lite\Cfg\Debug_Logs, especially

- install.log, install.move.log
- BadTrace.txt (if present)
- *.dmp files (if present) – include only those from the day when the issue occurred + from a few days before and after (these files can be quite large – if it is not possible to send them via email, you can use, e.g., uschozna.cz)
- *.log files from the PgmLog subfolder (since version 8.6)

4. Methods (*.met) used during the analysis, sequences (*.seq), measured chromatograms (*.prm), report styles (*.sty) and any other relevant files (e.g., exported data)

Clarity Lite function "Send Report by E-mail":

If the computer with Clarity Lite is connected to the internet, it is easier to let Clarity Lite collect the necessary information automatically. Note that it is required to have an e-mail client installed (other clients than Outlook may not work properly).

- Using the *Help – Send Report by E-mail* option from the main Clarity Lite window, a new e-mail message in your default e-mail client will be created, and the necessary common files will be automatically enclosed as an attachment.

If the size of the attachment is too large, consider sending the files using, e.g., uschovna.cz.

- Subsequently, you should fill in the description of the problem, enclose the relevant chromatogram (*.prm), sequence (*.seq), or report style (*.sty) files, and send the e-mail to the [DataApex technical support](#).