



## MARKES BENCHTOF

Clarity Control Module

ENG

Code/Rev.: M185/80A  
Date: 8/24/2020

Phone: +420 251 013 400  
Fax: +420 251 013 401  
clarity@dataapex.com  
www.dataapex.com

DataApex Ltd.  
Petrzilkova 2583/13  
158 00 Prague 5  
The Czech Republic

Clarity<sup>®</sup>, DataApex<sup>®</sup> and ▲<sup>®</sup> are trademarks of DataApex Ltd. Microsoft<sup>®</sup> and Windows<sup>™</sup> are trademarks of Microsoft Corporation.

*DataApex reserves the right to make changes to manuals without prior notice. Updated manuals can be downloaded from [www.dataapex.com](http://www.dataapex.com).*

Author: MP

# Contents

<b>1 Markes BenchTOF Control Module</b> .....	<b>1</b>
<b>2 Requirements</b> .....	<b>2</b>
2.1 Software requirements .....	2
2.2 Hardware requirements .....	2
<b>3 Installation Procedure</b> .....	<b>3</b>
3.1 Hardware - Wiring .....	3
3.2 Clarity Configuration .....	3
<b>4 Using the control module</b> .....	<b>5</b>
4.1 Method Setup - MS .....	5
4.2 Device Monitor .....	7
4.3 BenchTOF-dx Setup .....	9
<b>5 Report Setup</b> .....	<b>10</b>
<b>5 Troubleshooting</b> .....	<b>11</b>
5.1 Specific Problems .....	11

To facilitate the orientation in the **Markes BenchTOF** manual and **Clarity** chromatography station, different fonts are used throughout the manual. Meanings of these fonts are:

**Instrument** (blue text) marks the name of the window to which the text refers.

*Open File* (italics) describes the commands and names of fields in **Clarity**, parameters that can be entered into them or a window or dialog name (when you already are in the topic describing the window).

WORK1 (capitals) indicates the name of the file and/or directory.

*ACTIVE* (capital italics) marks the state of the station or its part.

The bold text is sometimes also used for important parts of the text and the name of the **Clarity** station. Moreover, some sections are written in format other than normal text. These sections are formatted as follows:

---

**Note:** Notifies the reader of relevant information.

---

**Caution:** Warns the user of possibly dangerous or very important information.

---

**■ Marks the problem statement or trouble question.**

**Description:** Presents more detailed information on the problem, describes its causes, etc.

**Solution:** Marks the response to the question, presents a procedure how to remove it.

# 1 Markes BenchTOF Control Module

This manual describes the setting of the **Markes BenchTOF** detector. The control module enables direct control of the instrument over LAN.



*Fig 1: Markes BenchTOF detector*

Direct control means that the MS detector can be completely controlled from the **Clarity** environment. Instrument method controlling the analysis conditions will be saved in the measured chromatograms.

## 2 Requirements

### 2.1 Software requirements

- **Clarity** Installation GC Control module (p/n A23) and MS-TOF Extension (p/n A37).

**Markes BenchTOF** requires **Microsoft .NET version 4** for operation. This version (or later) is already installed on majority of PCs. Nonetheless you will be notified during installation if your PC is missing **Microsoft .NET**. For complete list of .NET requirements, see the **.NET Framework System Requirements** on Microsoft web page.

---

**Caution:** Before installing **Clarity**, make sure that your **Windows** are updated to the latest version.

### 2.2 Hardware requirements

- LAN cable (p/n SK08).

---

**Note:** Cables are not part of **Clarity** Control Module. It is strongly recommended to order required cables together with the Control Module.

# 3 Installation Procedure

## 3.1 Hardware - Wiring

The **Markes BenchTOF** detector is controlled by LAN communication. It uses standard LAN cable.

## 3.2 Clarity Configuration





Fig 2: How to add Markes BenchTOF

- Start the **Clarity** station by clicking on the  icon on the desktop.
- Invoke the **System Configuration** dialog accessible from the **Clarity** window using the *System - Configuration...* command.
- Press the **Add** button ① (see **Fig 2** on pg 3.) to invoke the **Available Control Modules** dialog.
- You can specify the searching filter ② to simplify the finding of the driver.
- Select the **Markes BenchTOF** and press the **Add** ③ button.

The **BenchTOF-dx Setup** dialog will appear (see **Fig 3** on pg 4.).





Fig 3: BenchTOF-dx Setup

- Select the correct *IP Address* and *Port* numbers, then click the *Test* button to ensure the detector communicates correctly.

**Note:** Other fields from this dialog are described later in the description of the **BenchTOF-dx Setup** dialog.

- Press the **OK** button.

The **BenchTOF-dx** will appear in the *Setup Control Modules* list ④ of the **System Configuration** dialog.

- Drag and drop the **BenchTOF-dx** icon from the *Setup Control Modules* list ④ on the left side of the **System Configuration** dialog to the desired *Instrument* ⑤ tab on the right side ⑥ (or use the  button ⑦ to do so).

**Note:** The configuration dialog of the **Markes BenchTOF** detector (**BenchTOF-dx Setup**) can be displayed any time by double-clicking on its icon or using the *Setup* button.



## 4 Using the control module

New **MS** tab is created in the **Method Setup** dialog. It can be accessed by using the *Method - MS Control...* command from the **Instrument** window.

### 4.1 Method Setup - MS

The **Method Setup - MS** tab for the **Markes BenchTOF** detector allows to set the acquisition parameters of the **BenchTOF-dx** detector.




Fig 4: Method Setup - MS

#### Filament

Sets the working conditions for the filament in the **Markes BenchTOF** detector ion source.

##### Measurement Voltage

Sets the voltage on the electrodes during the measurement. The available values range between 0 and 3.8 Volts, with a step of 0.01 V.

##### Filament Delay

Sets the time after which will the voltage set on the filament be lowered to the voltage set in the *Standby Voltage* field to prevent excessive filament wear. The delay is set in seconds and may use values between 0 and 4000.

#### Temperature

Section allowing to set the temperatures for the different zones of the **Markes BenchTOF** detector.

**Ion Source Heating**

Sets the temperature for the ion source.

**Transfer Line Heating**

Sets the temperature for the transfer line.

**Ranges**

Sets the range of the signals in the scans using the *MS Range* field and range of the TIC signal using the *TIC Range* field. Both fields allow values between 0 and 1000000000 counts.

**Scans & Scansets**

Sets the data rate and scan size.


**Data Rate (Scansets per sec.)**

Defines the sample rate in terms of number of scans obtained per second. The size of the scan is defined using the *Start of Scan* and *End of Scan* fields, with the maximum scan size of 0 to 1747.7 amu. The *Start of Scan* value must be always at least 0.1 amu smaller than *End of Scan* value.

**Mass Representation**

Sets how the measured spectra are processed to create final spectra representation. The description of the fields is available in the **Markes BenchTOF** detector documentation made by **Markes**.

## 4.2 Device Monitor

The window with the detector status can be invoked by the *Monitor - Device Monitor* command from the [Instrument](#) window or using the  icon. It displays the actual values for detector parameters and allows to switch the detector between *OPERATE* and *STANDBY* states.




Fig 5: Device Monitor

**Caution:** Please head the instructions by **Markes** company on working with the **Markes BenchTOF** detector to prevent possible damage done to the detector.

### Temperature

Displays the *Set* and *Actual* temperatures on Ion Source and in the Transfer Line. *Set* values are taken for the current method, the *Actual* display the values as read from the detector.

### Filament

Displays the *Actual Voltage* and *Actual Current* on the filament as provided by the **Markes BenchTOF** detector.

### HW Info...

Button opening the [Hardware Configuration](#) dialog displaying the status of the detector including it's serial number.




Fig 6: Hardware Configuration

### Wake from Standby/Switch to Standby

Button switching the detector to and from the *STANDBY* state. If the detector is in the *STANDBY* state, no analysis can be performed in

---

**Clarity** Clarity. The detector switches down to the *STANDBY* state automatically after long inactivity, the time for the switch defined in the [BenchTOF-dx Setup](#) dialog.

## 4.3 BenchTOF-dx Setup

The **BenchTOF-dx Setup** dialog serves for the correct setting of the communication between **Clarity** and the **Markes BenchTOF** detector.




Fig 7: BenchTOF-dx Setup

### Connection

Section allowing to set the communication with the **Markes BenchTOF** detector. To establish the communication successfully the correct IP address must be entered into the *Host or IP Address* field. *Port Remote Interface* and *Port Data Interface* fields can be changed too, but same as with the IP address, the change must be also performed in the **Markes BenchTOF** to match the settings. Default values for both fields are 6392 and 6391, respectively.

### Test

Serves for detecting the **Markes BenchTOF** detector. If the communication attempt is successful, the serial number of the detector is displayed in the field next to the button.

### Standby Settings

Sets the time after which the **Markes BenchTOF** switches to the *STANDBY* state automatically during it's inactivity.

### External Trigger

Sets the necessary level of the detector's Start-in and Stop-in Input.

## 5 Report Setup

All of the detector settings accessible on the [Method Setup - MS](#) tab for the **Markes BenchTOF** detector are reported. To do so, the *Instrument Control* parameter on the [Method](#) tab of the [Report Setup](#) dialog must be checked.





Fig 8: Report Setup

# 5 Troubleshooting

## 5.1 Specific Problems

**In case your installation has frozen, it has been most probably caused by the fact that .NET Framework (dotNetFx40\_Full\_x86\_x64.exe as displayed by the picture below) has failed to install.**



- Solution:*
- 1) Using Windows **Task Manager**, end the installation process.
  - 2) Turn **OFF** Windows automatic updates and leave it turned **OFF** for the time being of the installation.
  - 3) Start Clarity installation again.
  - 4) After Clarity has been successfully installed, turn Windows automatic updates back **ON**.