

Clarity Demo

Clarity DEMO version allows you to try data acquisition and all related procedures as it provides data emulation.

Essential Features



- Try out any supported instrument without connecting a physical device
- Create and work with simulated data for a realistic experience
- Evaluate prepared demo chromatograms to explore Clarity's powerful data processing features

How to get Clarity Demo version

You can download our Clarity Demo version on www.dataapex.com after you register there. Then you simply install it and everything is ready to try out Clarity capabilities.

If you want to try the evaluation of your typical data, please let us know at sales@dataapex.com. We can prepare DEMO-ready chromatograms from AIA files (*.cdf) or ASCII text files containing signal values of your chromatograms.

Comparison Table

	 Clarity™ DEMO	 Clarity™
1000+ instruments available	✓ simulated control	✓ real connection
Acquisition	✓ simulated data	✓ real data
Data Import	✗	✓
Processing	✓ demo data	✓ real data
HW Key	✗	✓
Extensions	✓ all available	✓ purchased only

Extensions in Clarity DEMO

Clarity DEMO opens with four default projects: **GC**, **LC**, **GPC** and **LC-PDA**.

In addition, there are projects for the other extensions available in Clarity: **EA (Elemental Analysis)**, **NGA (Natural Gas Analysis)**, **DHA (Detailed Hydrocarbon Analysis)**, and **MS (Mass Spectrometry)**.

To launch them together with the prepared configuration, you can use our **Launch Manager utility**. Here, you select the demo project you want to launch and try out.