

Clarity – List of Currently Controlled Instruments (version 10.2)

GCs (p/n A23)

Producer	Name	Interface	Status
Agilent	4890D, 5890 II, 5890A	RS232	Ready
Agilent	6820, 6850 II	RS232 or LAN	Ready
Agilent	6850, 6890 (6890A), 6890 Plus, 6890N, 7820 (7820A), 7890 (7890A), 7890B, 8850, 8860, 8890, Intuvo 9000	LAN - ICF	Ready
Agilent ^{/1}	6850, 6890, 6890 Plus, 6890N	RS232	Ready
Ample Technology Center	ATC-6900 GC	LAN	Ready
Apix ^{/2}	ChromPix2, TwinPix	LAN	Ready
Dani	GC1000	RS232 + Colibrick	Ready
Dani ^{/2,3}	Master GC	RS232 or LAN	Ready
Ellutia ^{/2,4}	Series 200	RS232 + Colibrick	Ready
Ellutia ^{/2,4}	Series 500	LAN + Colibrick	Ready
Full ^{/2}	GC 9720 Plus, GC 9790 Plus	LAN	Ready
GOW-MAC	Series 816	LAN	Ready
GOW-MAC	Series 836	LAN - ICF	Ready
Konik ^{/5}	6000A	LAN	Testing
Labio Analytics Group S.a r.l. ^{/2}	GC11-UV	USB	Ready
Netel (India) Ltd. ^{/2}	Analyte 2900A, Chrom Lite 3000A	RS232 + Colibrick	Testing
Netel (India) Ltd. ^{/5}	Chrom Plus 6000	LAN	Testing
Shanghai Sunny Hengping Scientific Instrument ^{/2}	GC1290	LAN	Testing
Shimadzu ^{/2,6}	Brevis GC-2050, Nexis GC-2030	LAN - ICF	Ready
Shimadzu	GC14C, GC17A	RS232 + Colibrick	Ready
Shimadzu	GC2010, GC2010 Plus, GC2010 Pro, GC2014, GC2014C, GC2014C APC/AFC	RS232	Ready
Young In Chromass ^{/2}	ChroZen GC, ChroZen Plus GC	LAN	Testing
Young In Chromass ^{/2}	YL6100	RS232	Ready
Young In Chromass ^{/2}	YL6500	LAN	Ready

HPLC Systems (p/n A24)

Producer	Name	Interface	Status
Advion	Avant HPLC, Avant UHPLC	LAN	Ready
Agilent	1100 LC, 1200 LC, 1260 Infinity, 1290 Infinity II (see list of particular components)	LAN - ICF	Ready
Agilent ^{/1}	1100 system	GPIB	Ready

Agilent	1120 LC (G4286A, G4287A, G4288A, G4289A, G4290A), 1220 LC (G4286B, G4286C, G4287B, G4287C, G4288B, G4289B, G4290B, G4291B, G4291C, G4292B, G4292C, G4293B), 1220 VL (G4288C, G4289C, G4290C, G4293C, G4294B), 1260 Infinity II/III	LAN - ICF	Ready
Alltesta ^{/2}	Analyzer	RS232	Testing
Analytical Technologies Limited ^{/5}	iUHPLC Series 3000Plus	RS232	Testing
Antec	ALEXYS	RS232	Ready
Antec ^{/7}	ROXY	RS232	Ready
Axcend ^{/2}	Focus LC	LAN - ICF	Testing
CoruiTech ^{/2,8}	Rainbow, RainbowC	LAN	Testing
CVC Technologies ^{/5}	Atima	LAN	Testing
Ecom ^{/2}	Prepsystem (YMC)	RS232 or LAN	Ready
Esensing ^{/5}	MassChrom 1299	LAN + RS232	Testing
Hanon Instruments ^{/5}	LC7000, LC7300	LAN + RS232	Ready
Haosi Bio ^{/2}	ARIES LC 2400	RS485	Testing
Hitachi ^{/9,10}	Chromaster, Chromaster Plus	USB	Ready
Hitachi ^{/9,10,11}	Primaide	USB	Ready
Ingos ^{/2}	LC 5000	LAN	Testing
Interlab ^{/10}	MAESTRO	USB	Ready
Knauer ^{/2}	Azura, BlueShadow, Smartline, WellChrom	RS232	Ready
PerkinElmer ^{/2}	Flexar	LAN + RS232	Ready
PerkinElmer	NexSAR	RS232	Ready
RIGOL ^{/2}	L-3000	RS232	Testing
Shimadzu ^{/11}	LC-10, LC-20 Prominence	RS232 or LAN + corresp. controller (see HPLC Other)	Ready
Shimadzu ^{/2,3,6,11}	Nexera LC-40	LAN + corresp. controller (see HPLC Other)	Ready
Spark	LC1299, LC1299 Plus	RS232, USB or LAN	Testing
Sykam ^{/2}	S155, S500, S600	RS232	Ready
Watrex ^{/2}	Streamline	RS232 + USB	Ready
Welch	WiSys 5000	RS232 or LAN	Ready
Young In Chromass ^{/2}	ChroZen	LAN	Ready
Young In Chromass ^{/2}	YL9100, YL9300	LAN + RS232	Ready
Yusea-MS ^{/2}	LC2000	RS485 or LAN	Testing

HPLC Pumps (p/n A24)

Producer	Name	Interface	Status
Advion	Avant Pump	LAN	Testing
Agilent ^{/1}	1100 Binary, 1100 Isocratic, 1100 Quaternary	GPIB	Ready

Agilent	Binary (G1312A, G1312B, G1312C, G1376A, G2226A, G4220A, G4220B, G7112B, G7120A, G7132A, G7161A, G7161B), Isocratic (G1310A, G1310B, G1361A, G7110B), Quaternary (G1311A, G1311B, G1311C, G4204A, G5611A, G5654A, G7104A, G7104C, G7111A, G7111B, G7131A, G7131C) (all including selected pump valve clusters)	LAN - ICF	Ready
Agilent	Binary (G4302A, G4782A) (all including selected pump valve clusters)	LAN - ICF	Testing
Alltech/Grace	301, 426, 626, 726	RS232 - LC UNI	Ready
Analytical Technologies Limited	ALC2010, ALC2100, COLUMN-WASH	RS232 - LC UNI	Ready
Analytical Technologies Limited ^{/5}	P-3210U, PB-3220U, PQ-3245U	RS232	Testing
ANALYTIK Jena ^{/5}	PQ LC HPLC Pump 1125, PQ LC HPLC Pump 1130, PQ LC HPLC Pump 1132	RS232	Ready
Antec	LC100	RS232	Ready
Antec ^{/5}	LC110	RS232	Ready
Antec ^{/5}	LC110S	RS232 or LAN	Ready
ASI	500, 501, 520	RS232 - UNI Ruby	Ready
Beckman/Altex	110A, 110B, 114	Zebrick	Ready
Bio-Rad	1350	Zebrick	Ready
Bischoff	2250	RS232	Ready
Bischoff	3350	Zebrick	Ready
Chrom Tech	ISO-100, ISO-1000, ISO-1500, ISO-2000	RS232 - LC UNI	Ready
Chrom Tech	LD-Class	RS232 - UNI Ruby	Ready
Chromophor	Chromopeak 2010	RS232 - UNI Ruby	Ready
Chromsystems	CLC 300	Zebrick	Ready
Chromsystems ^{/5}	CLC 320, CLC 340	RS232	Ready
Chromsystems ^{/12}	CLC 330	RS232	Ready
CoruiTech ^{/2,8}	Corui Module BSM, Corui Module QSM	LAN	Testing
CQS ^{/5}	Climax 150 Pump, Climax 155 Pump, Climax S 1130	RS232	Ready
CXTH ^{/2}	P3000, P6000	RS232 or LAN	Ready
Ecom ^{/2}	Alpha 05, Alpha 10, Alpha 100, Alpha 250, Alpha 50, Beta 10, Beta 50, IOTA-250, IOTA-300, IOTA-400, IOTA-S-10, IOTA-S-50, IOTA-S-100, IOTA-S-200, IOTA-S-300, KAPPA10PP, LCP 4050, LCP 4100	RS232	Ready
Ecom	Alpha 10, Alpha 100, Alpha 50, Beta 10, Beta 50, KAPPA10, LCP 4100	RS232 - LC UNI	Ready
Ecom ^{/2}	ECP2002, ECP2010, ECP2011, ECP2050, ECP2100, ECP2300, ECP201L	RS232 or LAN	Ready
Ecom	LCP 4000, LCP 4020	Zebrick	Ready
Esensing	SPH1299	RS232, USB or LAN	Testing
Flom	011, 021, 022 (analog input optional)	Zebrick	Ready
Flom	301M, AI-12 series, UI-22 (UI series), UI-32	RS232 - UNI Ruby	Ready
Flux - Thermo ^{/2}	Allegro, Rheos 2000, Rheos 2200	RS232	Ready
Futecs ^{/2}	P-6000	RS232 - UNI Ruby	Ready
Gilson	302, 303, 305, 306, 307	GSIOC	Ready
Gynkotek	M300CS, M480	Zebrick	Ready
Hanon Instruments ^{/5}	LC7010, LC7011, LC7012	LAN	Ready

Hitachi ^{/9,10}	Chromaster 5110, 5160	USB	Ready
Hitachi	L-6000	via L-6200	Ready
Hitachi	L-6200, L-6210	RS232	Ready
Hitachi ^{/9,10,11}	Primaide 1110	USB	Ready
JASCO	PU-880	Zebrick	Ready
Knauer ^{/2}	10P/20P, 40P, 80P, P2.1S/P4.1S, P4.2S, PLATINblue P-1, S100, Smartline 1050	RS232 or LAN	Ready
Knauer ^{/2}	K-1000, K-1001, K-120, K-1800, Prep Pump 1800, K-501, Smartline S1000	RS232	Ready
Knauer	K-500, K-501, M-64	Zebrick	Ready
Knauer ^{/2}	P2.1L, P6.1L, P8.1L	LAN	Ready
Knauer	Smartline S1000	RS232 - LC UNI	Ready
Konik ^{/5}	560 Gradient Pump, 560 Isocratic Pump	RS232	Ready
Konik ^{/5}	580 Pump	RS232 or LAN	Ready
Kontron	320, 322, 325	RS232 - LC UNI	Testing
Kontron	420, 422, 525	RS232 - LC UNI	Ready
LCC Engineering	R 200, R 300, R 400	Zebrick	Ready
LKB - Amersham	LKB 3500, P-50, P-500	Zebrick	Ready
PerkinElmer ^{/2}	Flexar Binary LC Pump, Flexar Quaternary LC Pump	RS232	Ready
PerkinElmer	NexSAR pump	RS232 - UNI Ruby	Ready
Postnova Analytics	PN1130	Zebrick	Ready
Recipe	GP3000, IP3000	RS232	Ready
RIGOL ^{/2}	L-3200	RS232	Testing
Schambeck ^{/5}	S 9425	RS232	Ready
Schambeck	SDS 9404, SDS 9414 I	Zebrick	Ready
Schambeck	SFD 2100	RS232	Ready
Separlab ^{/2}	Separtrix PP03	RS232	Ready
Shimadzu	LC-10AD, LC-10ADvp, LC-10Ai, LC-10AS, LC-10AT, LC-10ATvp, LC-20AB, LC-20AD, LC-20ADXR, LC-20AT, LC-6AD, LC-8A	Shimadzu controller (see HPLC Other)	Ready
Shimadzu	LC-10ADvp, LC-10ATvp	OPT-USB	Ready
Shimadzu	LC-20AB, LC-20AD, LC-20ADXR, LC-20AT (in LC-10 emulation mode only)	OPT-USB	Ready
Shimadzu ^{/2,3,6}	LC-40B X3, LC-40B XR, LC-40D, LC-40D X3, LC-40D XR, LC-40D XS	Shimadzu controller (see HPLC Other)	Ready
Shimadzu	LC-6A	Zebrick	Ready
Shiseido Irica	IRICA	Zebrick	Ready
Spark	SPH1240	RS232, USB or LAN	Ready
Spark	SPH1299, SPH1299Q	RS232, USB or LAN	Testing
SSI	150, 1500, B2300, HF 300, LD-Class, LS-Class, LU-Class, M1-Class, Mighty-Mini, MX-Class, PR-Class, Prep 100, Prep 300, Q-Grad, Series I, Series II, Series III, UHP, VersaGrad, VersaGrad UHP	RS232 - UNI Ruby	Ready
SSI	222D	Zebrick	Ready
SunChrom ^{/5}	SunFlow 2010	RS232	Ready
Sykam	S1021, S1521, S1122	Zebrick	Ready
Sykam ^{/2}	S1125, S1130, S1132	RS232	Ready
Sykam	S2100	RS232	Ready
Techlab	Minipump1	RS232 - UNI Ruby	Ready
Teledyne	2300	Zebrick	Ready

Teledyne SSI	ReaXus LD-Class, ReaXus LS-Class, ReaXus LU-Class, ReaXus M1-Class, ReaXus MX-Class, ReaXus PR-Class	RS232 - UNI Ruby	Ready
TSP	Constametric III	Zebrick	Ready
Unicam	PU 4015	Zebrick	Ready
Unimicro ^{/2}	EasySep-3030 Pump	RS232	Testing
Varian	212-LC	RS232 - UNI Ruby	Ready
Waters	6000A, M 45, 501, 510, 515	Zebrick	Ready
Watrex	DeltaChrom P200, DeltaChrom SDS20, DeltaChrom SDS30	RS232 - LC UNI	Ready
Watrex	Mighty-Mini	RS232 - UNI Ruby	Ready
Watrex ^{/2}	Streamline P1m, Streamline P1s	USB	Ready
Welch ^{/5}	WelNova 6110 Isocratic pump, WelNova 6120 HPG pump, WelNova 6140 LPG pump	LAN	Ready
Young In Chromass ^{/2}	ChroZen Binary Pump, YL9110, YL9111, YL9112	LAN	Ready
Young In Chromass ^{/2}	ChroZen Pump	LAN	Testing
Young In Chromass	ChroZen UHPLC Pump	LAN	Testing
Yusea-MS ^{/2}	LC2000 Pump	RS485 or LAN	Testing
	Any pump with flow controlled by voltage or frequency	Zebrick	Ready
	Any pump with on-line flow controlled by text (ASCII) comm. protocol	RS232 - LC UNI	Ready

HPLC Detectors (p/n A24)

Producer	Name	Interface	Status
Agilent ^{/1}	1100 MWD, 1100 RID, 1100 VWD	GPIOB	Ready
Agilent	1200: ELSD	RS232	Ready
Agilent	ELSD (G4218A, G4260A, G4260B, G4261A, G4261B, G7102A), FLD (G1321A, G1321B, G1321C, G7121A, G7121B), MWD (G1365A, G1365B, G1365C, G1365D, G7165A), RID (G1362A, G7162A, G7162B), VWD (G1314A, G1314B, G1314C, G1314D, G1314E, G1314F, G7114A, G7114B)	LAN - ICF	Ready
Alltech/Grace	3300 ELSD	LAN	Ready
Analytical Technologies Limited ^{/5}	ASPD2600	RS232	Ready
Analytical Technologies Limited ^{/5}	UVD-3500U	RS232	Testing
ANALYTIK Jena ^{/2}	PlasmaQuant MS	LAN	Testing
Antec	DECADE Elite, DECADE Lite	USB or LAN	Ready
Antec	DECADE II, DECADE II MD, DECADE II SDC	RS232	Ready
Antec ^{/7}	ROXY Potentiostat	RS232	Ready
Chromsystems ^{/5}	CLC 420	RS232	Ready
CoruiTech ^{/2,8}	Corui Module UV Detector	LAN	Testing
CoruiTech ^{/2,8}	Rainbow detector	USB	Testing
CQS ^{/5}	Climax 150 IC, Climax 150 UV, Climax 155 IC, Climax S 3250 UV/Vis detector	RS232	Ready
CXTH ^{/2}	UV1000, UV3000, UV6000	RS232 or LAN	Ready

Ecom ^{/2}	Apatite, Flash 06 DAD, Flash 06 Single, Flash 10 DAD, Malachite, Opal, Sapphire, TOYFIX	RS232	Ready
Ecom ^{/2}	ECD2600, ECD2800	RS232 or LAN	Ready
Ecom ^{/2}	ECDM2000	RS232, USB or LAN	Ready
Ecom ^{/2}	LED20	RS232, USB or LAN	Testing
Ecom	RI2012	RS232 - UNI Ruby	Ready
Ecom	Topaz	RS232	Ready
Elysia-raytest ^{/2,13}	Gabi Nova	USB or LAN	Testing
Elysia-raytest ^{/2,13}	Gabi Star Radiodetector, Ramona Star	RS232	Ready
Elysia-raytest ^{/2,13}	POMO Radiodetector	USB - UNI Ruby	Ready
ESA ^{/14}	Coulochem III	USB	Ready
Esensing ^{/5}	ELSD	RS232	Testing
Esensing ^{/5}	MCH UV	LAN	Ready
Futechs ^{/2}	CD-6000, RI-6000	RS232 - UNI Ruby	Ready
Gilson	PrepELS™ II	RS232	Ready
Hanon Instruments ^{/5}	LC7020	LAN	Ready
Hanon Instruments ^{/5}	LC7070, LC7081	RS232	Ready
Hi Corporation	RI2012	RS232 - UNI Ruby	Ready
Hitachi ^{/9,10}	Chromaster 5410, 5420, 5440, 5450	USB	Ready
Hitachi ^{/9,10,11}	Primaide 1410	USB	Ready
IDEX (formerly ERC)	RefractoMax 520 Series	RS232 - UNI Ruby	Ready
Interlab	MAESTRO ELSD	RS232	Ready
Interlab	Omega FLD, Omega UV-VIS	USB	Ready
Japan Analytical Industry Co., LTD. ^{/2}	RI-700II	RS232	Ready
JASCO	FP-4020, FP-4025	RS232	Ready
Knauer ^{/2}	10D, 40D, 50D, CM2.1S, PLATINblue MW-1, S2520 VWD, S2550 VWD, UVD2.1L, UVD2.1S	RS232 or LAN	Ready
Knauer ^{/2}	Conductivity Monitor 2900, K-2600, S200 UV, S2300 RID, S2400 RID, Smartline S2500 VWD, WellChrom RI (K-2300, K-2301, K-2400, K-2401), WellChrom UV (K-2000, K-2001, K-2500, K-2501)	RS232	Ready
Knauer ^{/2}	Fluorescence Detector FLD2.1L	USB	Ready
Knauer ^{/2}	MWD2.1L, RID2.1L, VWD 2.1L	LAN	Ready
Knauer	UAD 6.1L, UAD 6.1L Prep	RS232 or USB	Ready
Konik ^{/5}	560 UV-Vis, 560/580 DAD	RS232	Ready
Konik	580 ELSD	RS232	Ready
Konik	580 RI	RS232 - UNI Ruby	Ready
Konik ^{/5}	580 UV-Vis	RS232 or LAN	Ready
Kontron	432	RS232 + Colibrick	Ready
Marion Research	Fixed Wavelength detector (255, 280, 415 nm)	USB	Ready
PerkinElmer ^{/2}	Flexar FX-UV/Vis, Flexar UV/Vis	RS232	Ready
PerkinElmer ^{/2,11}	NexION 1000, NexION 1100, NexION 2000, NexION 2200, NexION 300, NexION 350, NexION 5000	LAN	Ready
Recipe	EC6000	RS232 - UNI Ruby	Ready
Recipe ^{/5}	FW3000, UV3000	RS232	Ready
RIGOL ^{/2}	L-3500	RS232	Testing
RIGOL	L-3530	RS232	Ready

Runge ^{/2}	mikron 31, mikron 81	USB	Ready
Schambeck	RI 2000, RI 2012, S 2020, S 2030	RS232 - UNI Ruby	Ready
Schambeck ^{/5}	S 4245, ZAM 3000, ZAM 4000	RS232	Ready
Schambeck	S3210, S3240	RS232	Testing
Sedere	Sedex 100LT, Sedex 90LT, Sedex FP, Sedex LC	RS232 or USB	Ready
Sedere ^{/15}	Sedex 80LT	RS232	Ready
Sedere	Sedex 85LT	RS232	Ready
Semba Biosciences ^{/5}	Octave 4X, Octave 4XE	RS232	Ready
Shimadzu ^{/2,3,6}	CDD-10Avp, RF-20A, RF-20Axs, RID-20A, SPD-40, SPD-40V	Shimadzu controller (see HPLC Other)	Ready
Shimadzu	RF-10A, RF-10Axl, RF-20A, RF-20Axs, RID-10A, SPD-10A, SPD-10Ai, SPD-10AV, SPD-10AVi, SPD-10Avp, SPD-10AVvp, SPD-20A, SPD-20AV	Shimadzu controller (see HPLC Other)	Ready
Shimadzu ^{/5}	RF-10Axl, RF-20A, RF-20Axs	RS232	Ready
Shodex	RI-101, RI-102, RI-104	RS232 - UNI Ruby	Ready
Shodex	RI-501, RI-501EX, RI-502, RI-504	LAN - UNI Ruby	Ready
Softa	1300, 1400, 200S, 300S, 400	RS232	Ready
SunChrom	RefractoFlow 2030	RS232	Ready
SunChrom ^{/5}	SpectraFlow 2020	RS232 or LAN	Ready
Sykam	RI2000, S3585, S3590	RS232 - UNI Ruby	Ready
Sykam ^{/2}	S150, S3120, S3245, S3250	RS232	Ready
Sykam	S3210, S3240	RS232	Testing
Unimicro ^{/2}	EasySep-3030 UV Detector, EasySep-3030 UV/Vis Detector, UNIEX-7700 ELSD	RS232	Testing
VICI - Valco ^{/2}	TCD 3	RS232	Ready
VWR	VWR ELSD	RS232	Ready
Watrex ^{/5}	Streamline IC1	USB	Ready
Watrex ^{/5}	Streamline UV1	RS232	Ready
Wayeal ^{/2,8}	FLD3200	RS232 or LAN	Testing
Young In Chromass ^{/4,14}	ChroZen HPLC FLD, YL9140 FLD	USB	Ready
Young In Chromass ^{/2}	ChroZen HPLC RID-E, YL9170, YL9180	RS232	Ready
Young In Chromass ^{/2}	ChroZen HPLC-A UVD, ChroZen UHPLC UVD, YL9120	LAN	Ready
Young In Chromass	ChroZen RID, ChroZen UHPLC RID, YL9172 RID	LAN - UNI Ruby	Ready
Young In Chromass ^{/2}	ChroZen UVD	LAN	Testing
Young In Chromass	YL9181	RS232	Ready
Yusea-MS ^{/2}	LC2000 Detector	RS485 or LAN	Testing

PDA Detectors (p/n A24 + p/n A29)

Producer	Name	Interface	Status
Advion ^{/5}	Avant PDA Detector	USB or LAN	Ready
Agilent ^{/1}	1100: DAD, FLD	GPIB	Ready
Agilent	FLD (G1321A, G1321B, G1321C, G7121A, G7121B), G1315A, G1315B, G1315C, G1315D, G4212A, G4212B, G7115A, G7117A, G7117B, G7117C	LAN - ICF	Ready

Analytical Technologies Limited ^{/5}	DAD-3520U	RS232	Testing
ANALYTIK Jena ^{/5}	PQ LC PDA Detector	LAN	Testing
Axcend ^{/5}	FocusArray	LAN	Ready
Chromophor ^{/5}	PDA Detector	LAN	Testing
CoruiTech ^{/2,8}	Corui Module PDA	LAN	Testing
CQS ^{/5}	Climax S 3345 PDA Detector, Climax S 3350 PDA Detector	LAN	Testing
Ecom ^{/2}	BABY18DAD, TOY14DAD, TOY18DAD	RS232 or LAN	Ready
Ecom ^{/2}	ECDA2000, ECDA2800	USB or LAN	Ready
Ecom ^{/2}	Flash 12 DAD, FLASH 14 DAD, TOYDAD 400, TOYDAD 400 L, TOYDAD 600, TOYDAD 600 L, TOYDAD 800, TOYDAD 800 L	RS232	Ready
Esensing ^{/5}	MCH PDA	LAN	Ready
Hanon Instruments ^{/5}	LC7060	LAN	Ready
Hitachi ^{/9,10}	Chromaster 5430	USB	Ready
Hitachi ^{/9,10,11}	Primaide 1430	USB	Ready
Knauer ^{/2}	DAD2.1L, DAD6.1L, K-2700, K-2800, PLATINblue PDA -1, S2800 PDA, S2850 PDA	LAN	Ready
Knauer ^{/2}	Smartline S2600 MWD	RS232	Ready
Konik ^{/5}	580 PDA	USB or LAN	Ready
PerkinElmer ^{/2}	Flexar Plus PDA	LAN	Ready
Prince ^{/2}	Next DAD	RS232	Testing
RIGOL ^{/2}	L-3520	USB or LAN	Testing
Schambeck ^{/5}	S 4345, S 4350	LAN	Testing
Schambeck	S3210, S3240	RS232	Testing
Shimadzu	SPD-M10Avp, SPD-M20A, SPD-M30A	Shimadzu controller (see HPLC Other)	Ready
Shimadzu ^{/2,3,6}	SPD-M30A, SPD-M40A	Shimadzu controller (see HPLC Other)	Ready
SunChrom ^{/5}	SpectraFlow 2050	USB or LAN	Ready
Sykam	S3210, S3240	RS232	Testing
Sykam ^{/2}	S3345, S3350	LAN	Testing
Watrex ^{/5}	Streamline PDA	LAN	Testing
Young In Chromass ^{/2}	ChroZen HPLC-A PDA, ChroZen PDA, ChroZen UHPLC PDA, YL9160	LAN	Ready

LC-MS Detectors (p/n A38 + p/n A24)

Producer	Name	Interface	Status
Advion	Expression CMS	USB	Ready
Gilson	Verity 1920	USB	Ready
Konik	Q4 LCMS	USB	Ready
Microsaic	4000 MiD, 4500 MiD	LAN	Ready
QuadroCore ^{/2}	MCore	LAN	Testing
Young In Chromass	YL9900	USB	Ready

GC-MS Detectors (p/n A38 + p/n A23)

Producer	Name	Interface	Status
Konik ^{/5}	MSQ1	LAN	Testing
Shanghai Sunny Hengping Scientific Instrument ^{/2}	MSQ8100	LAN	Testing
Young In Chromass	ChroZen GC/MS, YL6900	LAN	Testing

GC-MS TOF Detectors (p/n A37 + p/n A23)

Producer	Name	Interface	Status
Markes International	BenchTOF	LAN	Ready

HPLC Other (p/n A24)

Producer	Name	Interface	Status
Agilent	G7180A	LAN - ICF	Ready
Bronkhorst ^{/5}	mini Cori-Flow	RS232	Ready
Ecom ^{/2}	ECLM2000	RS232, USB or LAN	Testing
Ecom	Panda-30 HID	USB	Ready
Knauer ^{/2}	ASM2.1L, ASM2.2L, IFU 2.1 LAN	LAN	Ready
Knauer ^{/2}	HPLC-Box, IF2, Manager 40M, Smartline Manager 5000	RS232	Ready
Knauer ^{/2}	IFU 2.1	USB	Ready
RotaChrom ^{/2}	rCPC	LAN	Testing
Shimadzu	CBM-20A, CBM-20ALite (for LC-10 or LC-20 system)	RS232 or LAN	Ready
Shimadzu ^{/2,3,6}	CBM-40, CBM-40 Lite, SCL-40 (for Nexera system)	LAN	Ready
Shimadzu	SCL-10Avp (for LC-10 system)	RS232	Ready
SISw	SSI pumps USB interface	USB - UNI Ruby	Ready
SISw ^{/2}	Systec Degasser Watchdog, USB pH monitor	USB - UNI Ruby	Ready
Testa ^{/2}	Flowmeter	USB - ICF	Testing
Unimicro ^{/2}	HVPS	RS232	Testing

GC Accessories (p/n A23)

Producer	Name	Interface	Status
Faculty of Science, Charles University in Prague	GCxGC Modulator	RS232	Ready

Autosamplers (p/n A26)

Producer	Name	Interface	Status
Advion	Avant Autosampler	LAN	Ready
Agilent ^{/1}	1100 ALS, 1100 TAS	GPIB	Ready
Agilent ^{/1}	6850 ALS	through GC6850	Ready

Agilent	7650, 7683A, 7683B, 7693A, 7697A, G1327A, G1368A, G1378A, G4226A, G5667A, G5668A, G7129A, G7129B, G7129C, G7137A, G7157A	LAN - ICF	Ready
Agilent ^{/1}	7673, 7673 II	RS232	Ready
Agilent	7673B, 7673C, 8697, 8697 XL, ALS 6850, G2258A, G3167A, G3167B, G4227A, G4303A, G4767A, G7158B, G7169B	LAN - ICF	Testing
Agilent ^{/1}	7683	through GC6890	Ready
Agilent	G1313A, G1329A, G1329B, G1367A, G1367B, G1367C, G1367D, G1367E, G1377A, G1389A (all including thermostat (G1330A, G1330B))	LAN - ICF	Ready
Agilent ^{/12}	G1888	LAN - ICF	Testing
Agilent	G2260A (all including thermostat (G1330A, G1330B))	LAN - ICF	Testing
Agilent	G7167A, G7167B, G7167C (all including thermostat (G7167-60201))	LAN - ICF	Testing
Agilent ^{/5}	PAL3 LSI 85 (G7366A), PAL3 RSI 120 (G7368A), PAL3 RSI 85 (G7367A), PAL3 RTC 120 (G7370A)	LAN - ICF	Testing
Alltech/Grace	570, 580	RS232	Ready
Ample Technology Center	ATC3000A, ATC3100A, ATC3200A	LAN	Ready
Analytical Technologies Limited ^{/5}	AS-3320U	RS232	Testing
Analytical Technologies Limited	AS2000	RS232	Ready
ANALYTIK Jena ^{/5}	PQ LC 5250, PQ LC 5300	RS232	Ready
Antec	AS 100	RS232	Ready
Antec ^{/5}	AS 110 Micro	RS232	Ready
Cetac	ASX-8000	USB	Ready
Chromsystems	CLC 200	RS232	Ready
Chromsystems ^{/5}	CLC 220	RS232	Ready
CoruiTech ^{/2,8}	Corui Module SamplerC, Corui Module SamplerU	LAN	Testing
CQS ^{/5}	Climax 150 AS, Climax S 5250, Climax S 5300	RS232	Ready
CTC Analytics	A200S	RS232	Ready
CTC Analytics	PAL (GC, LC, Combi), PAL-xt (LC-xt, HTC-xt, HTS-xt, HTX-xt, GC-xt, Combi-xt, Combi-xt extended)	RS232 or LAN	Ready
CTC Analytics	PAL LSI, PAL RSI, PAL RTC, PAL3 System	LAN - ICF	Testing
Dani	ALS1000	RS232	Ready
Dani ^{/2,3}	Master AS, Master DHS	RS232	Ready
Dani ^{/2,3}	Master SHS	RS232, USB or LAN	Ready
Dionamix ^{/2}	L-3320	RS232, USB or LAN	Testing
Ecom	AS96	RS232, USB or LAN	Ready
Eicom	AS-700	RS232, USB or LAN	Ready
ESA	542	RS232	Ready
Esensing	MCH Alias	RS232, USB or LAN	Ready
EST Analytical	AS120, Cobra	RS232	Ready
Hanon Instruments ^{/5}	LC7050	RS232	Ready
Hitachi ^{/9,10}	Chromaster 5210, Chromaster 5260, Chromaster 5280	USB	Ready
Hitachi ^{/9,10,11}	Primaide 1210	USB	Ready
HTA ^{/2}	2000H, 2000HT, 2100H, 2800T, 3000A, 3100A, 3200A	LAN - ICF	Testing

	HTA ^{/2}	HT1500L	USB	Ready
	HTA	HT2000H, HT2000HT, HT2100H, HT2800T, HT3000A, HT3100A, HT3200A, HT4000L	LAN	Ready
	HTA	HT200H, HT250D, HT280T (without SPME), HT300A, HT310A, HT300LV, HT310LV, HT300L, HT310L, HT800L	RS232	Ready
	Knauer	AS 3800, AS 3900, Optimas	RS232	Ready
	Knauer ^{/2}	AS 3950, AS 6.1L, PLATINblue AS-1	RS232 or LAN	Ready
	Knauer ^{/2}	AS 8.1L, Liquid Handler LH8.1	LAN	Ready
	Konik	AS 580	RS232, USB or LAN	Ready
	Kontron	360, 460, 465, 560	RS232 - UNI Ruby	Ready
	PerkinElmer ^{/2}	Flexar AS	RS232	Ready
	PerkinElmer	NexSAR Inert HPLC Autosampler	RS232, USB or LAN	Ready
	Postnova Analytics	PN5300	RS232	Ready
	Recipe	AS3000, AS5000	RS232	Ready
	Schambeck ^{/5}	S 6250, S 6300	RS232	Ready
	Schambeck	S5200	RS232	Ready
	Shimadzu	AOC-20i, AOC-20s, AOC-20d	RS232	Ready
	Shimadzu ^{/2,6}	AOC-20i, AOC-20s, AOC-20d	RS232 - ICF	Ready
	Shimadzu	AOC-5000 Plus	RS232 or LAN	Ready
	Shimadzu	SIL-10ADvp, SIL-10AF, SIL-10AXL, SIL-20A, SIL-20AC, SIL-20ACHT, SIL-20ACXR, SIL-20AHT, SIL-20AXR, SIL-HTA, SIL-HTC	Shimadzu controller (see HPLC Other)	Ready
	Shimadzu ^{/2,3,6}	SIL-40, SIL-40 XR, SIL-40C, SIL-40C X3, SIL-40C XR, SIL-40C XS	Shimadzu controller (see HPLC Other)	Ready
	Spark	Alias, Alias UHPLC (including User Program)	RS232, USB or LAN	Ready
	Spark	Basic Marathon (Plus, XT), Endurance, Midas, Triathlon	RS232	Ready
	Spark	Integrity	RS232, USB or LAN	Ready
	Spark ^{/2}	SP870	LAN - ICF	Testing
	Sykam	S5200	RS232	Ready
	Sykam ^{/2}	S5250, S5300	RS232	Ready
	VICI - Valco	Multi-position valve used as AS	RS232 - UNI Ruby	Ready
	Watrex ^{/5}	Streamline AS1, Streamline AS2	RS232	Ready
	Welch	WelNova 6210	RS232, USB or LAN	Ready
	Young In Chromass ^{/5}	ChroZen 2310	RS232	Testing
New	Young In Chromass ^{/5}	ChroZen 2330	LAN	Testing
	Young In Chromass ^{/2}	ChroZen AS, YL9150 Plus	LAN	Testing
	Young In Chromass	ChroZen UHPLC ALS, ChroZen3000A, ChroZen3100A, ChroZen3200A, YL2000H, YL2000HT, YL2100H, YL3000A, YL3100A, YL3200A	LAN	Ready
	Young In Chromass ^{/2}	YL3050A	RS232	Testing
	Young In Chromass ^{/5,16}	YL9150	RS232	Ready
	Young In Chromass	YL9151	RS232	Ready
	Zhengzhou Annuo ^{/2}	AS-3016, AS-3117	RS232	Testing
New	Zhengzhou Annuo ^{/2}	FAS-3823	LAN	Testing
		Any autosampler with the Ready input and Inject output (Inject marker)	Colibrick	Ready

HPLC Autosampler Accessories (p/n A26)

Producer	Name	Interface	Status
R-Biopharm	RIDACREST	RS232, USB or LAN	Ready
Spark	ACE, HPD, Symbiosis	RS232, USB or LAN	Ready

HPLC Column Ovens (p/n A24)

Producer	Name	Interface	Status
Advion	Avant Column Oven	LAN	Ready
Agilent ^{/1}	1100 TCC	GPIB	Ready
Agilent	G1316A, G1316B, G1316C, G7116A, G7116B, G7130A	LAN - ICF	Ready
Analytical Technologies Limited ^{/5}	CO-3400U	RS232	Testing
ANALYTIK Jena ^{/5}	PQ LC 4115	RS232	Ready
Chromsystems ^{/5}	CLC 360	RS232	Ready
Cluzeau Info-Labo	Croco-CIL, Gecko-2000, Igloo-CIL	RS232 - UNI Ruby	Ready
CoruiTech ^{/2,8}	Corui Module CHManager, Corui Module ColumnOven	LAN	Testing
CQS ^{/5}	Climax S 4115, Climax S 4120	RS232	Ready
Ecom ^{/2}	ECO2080, ECO2099	RS232 or LAN	Ready
Ecom	LCO 102	RS232 - UNI Ruby	Ready
Esensing	Mistral	RS232	Ready
Futecs ^{/2}	CT-6000	RS232 - UNI Ruby	Ready
Hanon Instruments ^{/5}	LC7030	LAN	Ready
Hitachi ^{/9,10}	Chromaster 5310, Chromaster Ultra RS 6310	USB	Ready
Hitachi ^{/9,10,11}	Primaide 1310	USB	Ready
Knauer ^{/2}	CT2.1, PLATINblue T-1, S4000, S4050	RS232 or LAN	Ready
Knauer ^{/2}	JetStream	RS232	Ready
Konik ^{/5}	560 Oven	RS232	Ready
Konik ^{/5}	580 Oven	RS232 or LAN	Ready
PerkinElmer ^{/2}	Flexar LC Column Oven	RS232	Ready
PerkinElmer	NexSAR Column Oven	RS232	Ready
R-Biopharm	RIDACREST Column Oven	RS232	Ready
Recipe	CH3000	RS232 - UNI Ruby	Ready
Recipe	HT4000	RS232 or LAN	Ready
RIGOL ^{/2}	L-3400	RS232	Testing
Schambeck ^{/5}	S 5120	RS232	Ready
Shimadzu	CTO-10A, CTO-10AC, CTO-10ACvp, CTO-10ASvp, CTO-10Avp, CTO-20A, CTO-20AC	Shimadzu controller (see HPLC Other)	Ready
Shimadzu ^{/2,3,6}	CTO-30A, CTO-40C, CTO-40S	Shimadzu controller (see HPLC Other)	Ready
Spark	Mistral	RS232 or LAN	Ready
Spark ^{/2}	SP240	LAN - ICF	Testing
SSI	605 COC	RS232 - UNI Ruby	Ready
SunChrom ^{/5}	SunTherm 2070	RS232 or LAN	Ready

Sykam ^{/2}	S4115, S4120	RS232	Ready
Techlab	K-7	RS232 - UNI Ruby	Ready
Thermotechnic Products ^{/5}	JetStream 2 Plus	RS232	Ready
Unimicro ^{/2}	EasySep-3030 Thermostat, Unimicro Thermostat	RS232	Testing
VICI - Valco ^{/2}	Fast Temperature Programmer	RS232	Ready
Watrex	CTC 100	RS232 - UNI Ruby	Ready
Watrex ^{/2}	Streamline CT1	USB	Ready
Young In Chromass ^{/2}	ChroZen Column Oven	LAN	Testing
Young In Chromass ^{/2}	ChroZen HPLC-A CC, ChroZen UHPLC CC, YL9130	LAN	Ready
Yusea-MS ^{/2}	LC2000 Thermostat	RS485 or LAN	Testing

Fraction collectors (p/n A24)

Producer	Name	Interface	Status
Advantec	CHF122SC	RS232	Ready
Agilent	G1364A, G1364B, G1364C, G1364D, G1364E, G1364F, G5664A, G5664B, G7159B, G7166A	LAN - ICF	Testing
Bio-Rad	BioFrac	Colibrick + FC GP	Ready
Büchi ^{/5,17}	C-660	RS232	Ready
Ecom ^{/2}	ECF2096	RS232 or LAN	Ready
Ecom	ECV2010	RS232	Ready
Ecom ^{/2}	FRC ECVA2000	RS232 or USB	Ready
Gilson	203B, 204	RS232	Ready
Japan Analytical Industry Co., LTD. ^{/2}	FC-7000	RS232	Ready
Knauer ^{/2,17}	FC 6.1	LAN	Ready
Knauer ^{/2,17}	S3000, S3050	RS232 or LAN	Ready
Labomatic Instruments ^{/5,17}	Labocol VARIO-4000, Labocol VARIO-4000 Plus	RS232 or LAN	Ready
LKB Amersham	Frac100	Colibrick + FC GP	Testing
Rheodyne	MX Series II valves used as FC	USB	Ready
Separlab ^{/2}	Separflow FC	RS232	Ready
Sykam	S6410	Colibrick + FC GP	Ready
Teledyne	FOXY 200, FOXY Jr., FOXY R1, FOXY R2	RS232	Ready
Teledyne ^{/5,17}	FOXY Jr.	RS232	Ready
Teledyne ^{/5,17}	FOXY R1, FOXY R2	RS232 or LAN	Ready
VICI - Valco	Multi-position valve used as FC	RS232 - UNI Ruby	Ready
Waters	Fraction Collector II	RS232	Ready
Watrex	DeltaChrom FC 500	Colibrick + FC GP	Ready
Watrex ^{/2}	StreamLine FC	USB	Ready
Young In Chromass	YL FC1, YL FC2	RS232	Ready
	Any fraction collector using the Next and Collect/Waste events	Colibrick - FC-GP	Ready

Valves (p/n A24 or A23)

Producer	Name	Interface	Status
Agilent	G1156A, G1157A, G1158A, G1158B, G1159A, G1160A, G1162A, G1163A, G1170A, G4243A, G7170B, G9322A	LAN - ICF	Testing
ANALYTIK Jena ^{/5}	PQ LC 6510	RS232	Ready
Chromsystems ^{/5}	CLC 230	RS232	Ready
Ecom ^{/2}	ECVA2000	RS232, USB or LAN	Ready
Knauer ^{/2}	K-6, K-12, K-16, V2, V6, V12, V16, V2.1S, VU4.1	RS232 or LAN	Ready
PerkinElmer	NexSAR Switching valve	RS232	Ready
Rheodyne	LabPro Series valves	RS232 - UNI Ruby	Ready
Rheodyne	MX Series II valves	USB	Ready
Rheodyne	Titan Series valves	RS232	Ready
Shimadzu ^{/2,3,6}	FCV-0206, FCV-0206H, FCV-0206H3, FCV-0607, FCV-0607H, FCV-0607H3, FCV-DR (all through selected Nexera column ovens only)	LAN - ICF	Ready
Shimadzu	FCV-10AL, FCV-10ALvp, FCV-11AL, FCV-11ALS, FCV-15AL, LPGE Unit (all through selected LC-10/20 pumps only)	Shimadzu controller (see HPLC Other)	Ready
Shimadzu ^{/2,3,6}	FCV-11AL, FCV-11ALS, LPGE Unit (all through selected Nexera pumps only)	LAN - ICF	Ready
Shimadzu	FCV-12AH, FCV-14AH (all through selected LC-20 column ovens only)	Shimadzu controller (see HPLC Other)	Ready
Sykam ^{/2}	S6510	RS232	Ready
Upchurch	V-2500	RS232	Ready
VICI - Valco	Any valve mounted on two-position microelectric actuator, multiposition microelectric actuator or universal actuator.	RS232	Ready

Amino Acid Analyzers

Producer	Name	Interface	Status
Ecom ^{/5}	ECA AMINO 500	LAN	Testing
Ingos ^{/2}	AAA 500	LAN	Testing
Sykam ^{/4}	S4300 Amino Acid Reaction Module System	RS232	Ready
Sykam ^{/2,4}	S633	RS232	Testing

Capillary Electrophoresis (p/n A24 + p/n A31)

Producer	Name	Interface	Status
CMP Scientific ^{/5}	ECE-001	LAN	Testing
Prince ^{/2}	Next CE	LAN	Testing

Balances for Clarity EA Extension (p/n A30)

Producer	Name	Interface	Status
Mettler Toledo	AT, MT, UMT, Excellence	RS232	Ready

Sartorius	CP, GC and GP series, Cubis series, M2P	RS232	Ready
-----------	---	-------	-------

Other (p/n A24 or A23)

Producer	Name	Interface	Status
DataApex	Column Usage Monitor	UNI Ruby	Ready

Analog GCs

Producer	Name	Interface	Status
ABB	PGC5000B, PGC5009	Colibrick	Ready
Agilent	4890A	Colibrick	Ready
Agilent	5880 (analog output optional)	Colibrick	Ready
AIT	FXI Series5	Colibrick	Ready
AMA Instruments	GC 5000	Colibrick	Ready
ASDevices	GCSense, miniMOv	Colibrick	Ready
CHEMTRON	2045	Colibrick	Ready
Chromatec	Kristall 2000M	Colibrick	Ready
Chromatotec	GC 866	Colibrick	Ready
Chrompack	CP9001, CP9002	Colibrick	Ready
DPS	Companion, DPS 600	Colibrick	Ready
Emerson	Rosemount 370XA, Rosemount 700XA	Colibrick	Ready
Global Analyser Solutions	CompactGC (analog output optional)	Colibrick	Ready
GOW-MAC	20 Series, 23-550-1 and 23-750-1, Series 100, Series 1300, Series 1400, Series 200BTEX, Series 200SA, Series 23-550-BG, Series 350, Series 400, 400-P, 400-GP, Series 400MGA, Series 50, Series 5200, Series 580, Series 590, 592, Series 5900, Series 600, Series 7100 & 7100CR, Series 8100, Series 900, Series AR710 HFADD	Colibrick	Ready
INFITEK	GC-900A, GC-906A	Colibrick	Ready
Inrag	I-Graph XPX	Colibrick	Ready
LABIO	GC 72	Colibrick	Ready
LDetek	LD8000, Multidetek2	Colibrick	Ready
LP	Chrom 4, Chrom 5, Chrom 61	Colibrick	Ready
Makarand	MAK 2191	Colibrick	Ready
NETEL	MICRO 9100, OMEGA QC	Colibrick	Ready
Nexus Instruments	M3	Colibrick	Ready
NTS	IATROSCAN MK-5, IATROSCAN MK-6	Colibrick	Ready
Nucon	5700, 5765, 5765 EPC	Colibrick	Ready
Orthodyne	Argon, Line 500	Colibrick	Ready
Peak Laboratories	PEAK PERFORMER	Colibrick	Ready
PerkinElmer	Autosystem, Autosystem XL, Clarus 400, Clarus 480, Clarus 500, Clarus 580, Clarus 590, Clarus 600, Clarus 680, Clarus 690, PE 8700, Sigma	Colibrick	Ready
PID Analyzers	322	Colibrick	Ready
SAKS	SEPmaster - 1000	Colibrick	Ready
Scion	SCION 436, SCION 450, SCION 456, SCION 8300, SCION 8500 (analog output optional)	Colibrick	Ready

Servomex	SERVOPRO DF-750, SERVOPRO NanoChrome ULTRA	Colibrick	Ready
Shimadzu	GC-14A, GC-14B, GC-2025, GC-8, GC-8A, GC-9A, Tracera	Colibrick	Ready
SRA	R990 (analog output optional)	Colibrick	Ready
Teledyne	Falcon GC (analog output optional)	Colibrick	Ready
Thermo Carlo Erba	6000 Vega, HRGC 8000, HRGC 8000 TOP, HRGC Mega 2, HRGC Mega 5300	Colibrick	Ready
THERMO Finnigan	GC 9001	Colibrick	Ready
Thermo Fisher	Model 42i, Trace 1110, TRACE 1300, TRACE 1310	Colibrick	Ready
Thermo Fisons	FOCUS, Trace GC 2000, Trace GC Ultra	Colibrick	Ready
Thermo India	GC 8610	Colibrick	Ready
Unicam	104 GC, PU 204, PU 4550, PU 4600	Colibrick	Ready
UnitGC	EmitChrom	Colibrick	Ready
Varian	3300, 3400, 3600, 3700, 3350, 3380, 3800, 3900, 450-GC	Colibrick	Ready
Varian	430-GC (analog output optional)	Colibrick	Ready
Yokogawa	GC8000	Colibrick	Ready

Analog HPLC Detectors

Producer	Name	Interface	Status
Adelis	Zetalif LEDIF, Zetalif LIF	Colibrick	Ready
Agilent	1046, 1047, 1049A, 1050 VWD	Colibrick	Ready
Alcott	795	Colibrick	Ready
Applied Biosystems	Kratos Spectroflow 757, Kratos Spectroflow 783, Kratos SpectroFlow 980	Colibrick	Ready
Aquilon	102M, 104.1M, 121M	Colibrick	Ready
ASDevices	Cubed	Colibrick	Ready
BASi	EC epsilon, LC epsilon, LC-3C, LC-3D, LC-3E, Petit Ampère, LC-4B, LC-4C	Colibrick	Ready
Beckman	Pace	Colibrick	Ready
Berthold	FlowStar LB 509, FlowStar LB 513, FlowStar LB 514, fLumo, HERM LB 500	Colibrick	Ready
Bio-Rad	1755	Colibrick	Ready
Bischoff	BioQuant PAM 2, LAMBDA 1010, RI 8120	Colibrick	Ready
Brechbühler	Sniffer 9000	Colibrick	Ready
Büchi	ELS Detector C-650	Colibrick	Ready
Carroll & Ramsey	Model 105S-1	Colibrick	Ready
China Yangzhou	WD-RI100	Colibrick	Ready
D-Star	DFW-20, DFW-20C, DFW-21C, DVW-10	Colibrick	Ready
Dionex	AD20, MWD-3000	Colibrick	Ready
Dionex	FLD-3100, FLD-3400RS, VWD-3100, VWD-3400RS (analog output optional)	Colibrick	Ready
Eckert & Ziegler	HPLC Scan	Colibrick	Ready
Ecom	LCD 2040, LCD 2070.1, LCD 2071.2, LCD 2071.3, LCD 2073A, LCD 2082, LCD 2083.61, LCD 2083.7, LCD 2084, LCD 2563	Colibrick	Ready
Eicom	ECD-700	Colibrick	Ready
ESA	Chromachem, Corona, Coulochem II	Colibrick	Ready
GBC	LC1200, LC1205, LC1210, LC1245, LC1260	Colibrick	Ready

Gilson	112, 151, 152, 155, 156, 157, 159	Colibrick	Ready
Grace Alltech	2000, 2000 ELSD, 500, 800	Colibrick	Ready
Hamamatsu	Photon's PMT	Colibrick	Ready
Hitachi	L-2400, L-2420, L-2480, L-2485, L-2490, Primaide 1440 (analog output and keyboard optional)	Colibrick	Ready
Hitachi	L-2450, L-2455 (single wavelength output only, analog output and keyboard optional)	Colibrick	Ready
Hitachi	L-3500A, L-4000, L-7400, L-7420, L-7480, L-7485, L- -7490	Colibrick	Ready
Hitachi	L-7450, L-7455 (single wavelength output only)	Colibrick	Ready
IBZ	CHIRALYSER-MP	Colibrick	Ready
IN/US Systems	β -RAM	Colibrick	Ready
Istech	TraceDec	Colibrick	Ready
JAI	RI-50s, UV-370 NEXT, UV-3702, UV-3740, UV-600 NEXT	Colibrick	Ready
JASCO	3177, CD-1595, CD-2095, CD-4095, CL-2027, FP -2020, FP-2025, FP-821, FP-920, OR-2090, RI-2031, UV-1575, UV-2070, UV-2077, UV-975	Colibrick	Ready
KNAUER	2000 ES, ClinLab EC3000, LB 509	Colibrick	Ready
Kontron	330, 332, 335, 430, 535, SFM-25	Colibrick	Ready
LabAlliance	201, 500, 504, 506	Colibrick	Ready
LDetek	Plasmatetek, PlasmaDetek 3	Colibrick	Ready
LKB Amersham	2138, UV-900, UV-M II	Colibrick	Ready
LP	RIDK100, RIDK101, RIDK102	Colibrick	Ready
Metrohm	Bioscan 817	Colibrick	Ready
OI Analytical	5383 PFPF	Colibrick	Ready
PAC	SeNSE	Colibrick	Ready
PerkinElmer	200 RID, 200 UV/VIS, 785A	Colibrick	Ready
Polymer Laboratories	ELS 1000, ELS 2100	Colibrick	Ready
Postnova	PN3150	Colibrick	Ready
Precision	IOTA 2	Colibrick	Ready
PSA	PSA Millenium Merlin	Colibrick	Ready
Quant	NQAD	Colibrick	Ready
Reach devices	RD4-250-280-RI	Colibrick	Ready
Recipe	EC3000	Colibrick	Ready
Repligen	Spectra/Chrom 280	Colibrick	Ready
SEDERE	60LT, 754, 75LT	Colibrick	Ready
Shimadzu	CDD-6A, RF-535, RF-551, RID-6A, SPD-6A	Colibrick	Ready
Shodex	CD-200, RI-201, RI-201H	Colibrick	Ready
Sievers	355 SCD	Colibrick	Ready
SoFTA	Model 2300, Model 300	Colibrick	Ready
SunChrom	SpectraFlow 501	Colibrick	Ready
Techcomp	ELSD180	Colibrick	Ready
Teledyne	UA-6	Colibrick	Ready
THERMO Finnigan	SpectraSystems RI-150	Colibrick	Ready
Thermo Fisher	FOCUS 2000, SM 4100	Colibrick	Ready
TSP	EC 2000, FL 2000, SM 3200, SP100, Spectromonitor D, UV 100, UV 1000, UV 2000, UV 3000	Colibrick	Ready
TSP/Thermo Fisher	FL 3000	Colibrick	Ready
Unicam	PU 4020, PU 4025, PU 4110	Colibrick	Ready

Varian	310, 356-LC, 380-LC, 385-LC, ProStar 325, ProStar 340, ProStar 345, ProStar 355, ProStar 363, ProStar 370	Colibrick	Ready
Visacon	VUV-14, VUV-14D, VUV-24, VUV-40F	Colibrick	Ready
Viscotek	VE3580	Colibrick	Ready
Waters	2410, 2414, 2420, 2424, 2465, 2475, 2487, 2489, 410, 431, 432, 464, 474, 480, 481, 484, 486, 490, 490E	Colibrick	Ready
Watrex	DeltaChrom UV200, DeltaChrom VD 400	Colibrick	Ready
Wyatt	Optilab rEX	Colibrick	Ready
Xylem	4430 PID	Colibrick	Ready
Yamazen	RI-31	Colibrick	Ready

Analog Amino Acid Analyzers

Producer	Name	Interface	Status
Biochrom	Biochrom 20, Biochrom 30+	Colibrick	Ready

Analog Elemental Analyzers

Producer	Name	Interface	Status
Costech/NC Technologies	8070 Air CO2, ECS-4010, ECS-4024, ECS-8020, ECS-8040	Colibrick	Ready
Thermo Carlo Erba	NA 1500, NC 2500	Colibrick	Ready

Other Analog Devices

Producer	Name	Interface	Status
Prince	PrinCE 500	Colibrick	Ready

Interface types:

RS232	=	Serial communication. Specify the type of the cable or instrument.
RS485	=	Serial communication. Specify the type of the cable or instrument.
LAN	=	Network communication. Specify the type of the cable or instrument.
USB	=	USB communication. Specify the type of the cable or instrument.
GPIB	=	hardware: Interface card for older Agilent 1100 system, Part No.: IAG11 (PCI interface) or IAG12 (USB interface)
GSIOC	=	hardware: Converter for Gilson serial communication, Part No.: IGLN1
Zebrick	=	hardware: External LC Control device (4x D/A & 4x D/F converter), Part No.: B30
OPT-USB	=	hardware: Optical to USB converter, Part No.: OPT01
LC UNI	=	software: Universal LC control module for online pumps, Part No.: A24
FC GP	=	software: General Purpose Fraction Collector control module for Clarity, Part No.: A24
UNI Ruby	=	software: Universal control module for various instruments, Part No.: A24
ICF	=	software: Agilent ICF control Libraries embedded in Clarity, Part No.: A24, A26, A29 depending on the type of controlled instruments
Colibrick	=	hardware: A/D converter also capable of sending and receiving TTL signals governed by Clarity. Part No.: U31, U32, U34

Status:

Ready	=	the module has been subject to extensive testing and has been officially RELEASED.
Testing	=	the module is already functional and its development is in stage of external testing often called public BETA. It is possible to get this control module and help us with fine tuning the driver for your specific configuration.

Notes:

1. *Control for compatibility with obsolete communication method.*
2. *Developed by instrument's producer.*
3. *Exclusively available from the instrument's producer – needs a release code / connection key.*
4. *Exclusively available from the instrument's producer in the particular OEM version.*
5. *Developed by third party.*
6. *Exclusively available in selected territories through companies approved by Shimadzu company.*
7. *Only models with Firmware 5.11 and lower.*
8. *Exclusively available from standalone installer.*
9. *Exclusively available in selected territories through companies approved by Hitachi company.*
10. *One module in the system must have Interface Card installed.*
11. *Only limited number of these systems on one Clarity station is possible - see appropriate manual or product page.*
12. *With limitations, contact support@dataapex.com*
13. *Selected models only - Check with Elysia-raytest company.*
14. *Limited to certain OS versions, see product page for details.*
15. *Only models with serial option - Check with Sedere company.*
16. *Only models with 840.100 firmware family (manufactured after 2010).*
17. *Available from Knauer company only (requires Knauer FRC option).*