

SW Protection - Clarity SW since version 2.4

The Clarity SW since version 2.4 (Clarity p/n C50, Clarity Lite p/n C40, Clarity Offline p/n C59) will be protected by a dongle (hardware key, either USB or LPT). Today's SW protection (CSW 17, CSW 32 and Clarity up to version 2.3) is based on combination of User Code and a specific A/D converter.

Both distributors and customers will benefit from higher flexibility. The essential difference and advantage of the dongle protection is ability to use any A/D converter* with any Clarity SW (since version 2.4.) There will be simplified ordering procedure without the necessity of specifying A/D cards for SW. Distributors will be able to keep stocks easier now.

* Please refer to the compatibility table reflecting various operating systems.



Rockey LPT and USB hardware keys (dongles)

DataApex will distribute USB dongles by default. LPT dongles will be available upon request. Every dongle includes information on S/N of Clarity SW and is accompanied by corresponding user code. IMPORTANT Please note that a loss of a dongle (hardware key) is considered as a loss of the license right to use Clarity SW.

A dongle installation will be a part of the Clarity SW installation (will be done automatically) and will be described in the *Getting Started* manual.

Download the Getting Started manual (M006) from www.dataapex.com section downloads

Clarity ver. 2.4 release date is February 1st, 2005. All Clarity users may update to Clarity version 2.4 for a fee (refer to Price list 2005).

Users of Clarity SW ordered after August 1st 2004 will receive this update free of charge (upon request).

For more information contact sales@dataapex.com